

# Goodman Acquisition Camera (GACAM) Cheat Sheet (Dec 2015)

**1)** Start GACAM software

**2)** From Goodman GUI, place slit mask IN

**3)** Center crosshair on slit. Select "Lock crosshair"

**4)** Withdraw slit mask (in Goodman GUI).  
Unlock Slit-Box un pop-up menu.

**5)**

**6)**

**7)** Unselect the "Show crosshair" option to see more clearly the target in the slit center.

**8)** Move the GACAM arm to the OUT position.  
Start your science exposure.

## Adjusting the exposure time and display scale:

Use the **Exposure** (left) and **Scale** (right) sub-windows to adjust the exposure time and minimum/maximum scaling values depending on the conditions of your target and sky. Also, experiment with the display scaling options in the second and third block of the pop-up menu.

Finally, you can always swap and use the slit box to define your slit center, and the crosshair for picking up a specific location on your target.